

BTB/POZ domain-containing protein 2 Polyclonal Antibody

Catalog No: #42422

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	BTB/POZ domain-containing protein 2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total BTB/POZ domain-containing protein 2 polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant human BTB/POZ domain-containing protein 2 protein
Target Name	BTB/POZ domain-containing protein 2
Other Names	BTBD2
Accession No.	Swiss-Prot#: Q9BX70
Uniprot	Q9BX70
GenEID	55643;
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:20 - 1:200

Images



Immunohistochemical analysis of paraffin-embedded human kidney using #42422 at dilution of 1:50.

Background

Interacts with topoisomerase 1 and with TRIM5 isoform Delta.

References

[1]"Characterization of BTBD1 and BTBD2, two similar BTB-domain-containing Kelch-like proteins that interact with topoisomerase I." Xu L., Yang L., Hashimoto K., Anderson M., Kohlhagen G., Pommier Y., D'Arpa P. BMC Genomics 3:1-1(2002) [PubMed: 1181802]

Note: This product is for in vitro research use only