

Activator of 90 kDa heat shock protein ATPase homolog 1 Polyclonal Antibody

Catalog No: #42435

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	Activator of 90 kDa heat shock protein ATPase homolog 1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	WB IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total Activator of 90 kDa heat shock protein ATPase homolog 1 polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant human Activator of 90 kDa heat shock protein ATPase homolog 1 protein
Target Name	Activator of 90 kDa heat shock protein ATPase homolog 1
Other Names	AHSA1
Accession No.	Swiss-Prot#: O95433
Uniprot	O95433
GeneID	10598;
Calculated MW	37kd
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

Application Details

Western blotting: □1:500 - 1:1000

Immunohistochemistry: 1:20 - 1:200

Images



All lanes : Activator of 90 kDa heat shock protein ATPase homolog 1 antibody at 2ug/ml

Lane 1 : EC109 whole cell lysate


Lane 2 : 293T whole cell lysate

Secondary

Goat polyclonal to Rabbit IgG at 1/15000 dilution

Predicted band size : 37 kDa

Observed band size: 45kDa



Immunohistochemical analysis of paraffin-embedded human spleen using #42435 at dilution of 1:100.

Background

Cochaperone that stimulates HSP90 ATPase activity By similarity. May affect a step in the endoplasmic reticulum to Golgi trafficking.

References

[1] "Gene expression profiling in the human hypothalamus-pituitary-adrenal axis and full-length cDNA cloning."Hu R.-M., Han Z.-G., Song H.-D., Peng Y.-D., Huang Q.-H., Ren S.-X., Gu Y.-J., Huang C.-H., Li Y.-B., Jiang C.-L., Fu G., Zhang Q.-H., Gu B.-W.,

Note: This product is for in vitro research use only