

C-X-C motif chemokine 13 Polyclonal Antibody

Catalog No: #42450

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	C-X-C motif chemokine 13 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total C-X-C motif chemokine 13 polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant Human C-X-C motif chemokine 13 protein
Target Name	C-X-C motif chemokine 13
Other Names	Angie, B cell-attracting chemokine 1, B lymphocyte chemoattractant, CXC chemokine BLC, Small-inducible cytokine B13, CXCL13, BCA1, BLC, SCYB13
Accession No.	Swiss-Prot#: O43927
Uniprot	O43927
GeneID	10563;
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:20 - 1:200

Images



Immunohistochemical analysis of paraffin-embedded human liver using #42450 at dilution of 1:100.

Background

Chemotactic for B-lymphocytes but not for T-lymphocytes, monocytes and neutrophils. Does not induce calcium release in B-lymphocytes.

References

[1]AB-cell-homing chemokine made in lymphoid follicles activates Burkitt lymphoma receptor-1.Gunn M.D., Ngo V.N., Ansel K.M., Ekland E.H., Cyster J.G., Williams L.T.Nature 391:799-803(1998) [2]B cell-attracting chemokine 1,a human CXC chemokine expressed

Note: This product is for in vitro research use only