

# Porcine Immunoglobulin A (IgA) ELISA Kit

Catalog No: #EK14659



Package Size: #EK14659-1 48T #EK14659-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Porcine Immunoglobulin A (IgA) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Sus scrofa; Porcine Pig
Storage	12 months

## Application Details

Detect Range:31.2-2000pg/mL

Sensitivity:13.2pg/mL

Sample Type:serum, plasma or other biological fluids.

Sample Volume: 100 µL

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only