

# Neurotrophin-4 Polyclonal Antibody

Catalog No: #42517

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	Neurotrophin-4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total Neurotrophin-4 polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant human Neurotrophin-4 protein
Target Name	Neurotrophin-4
Other Names	NT-4, Neurotrophin-5, NT-5, Neurotrophic factor 4, NTF4, NTF5
Accession No.	Swiss-Prot#: P34130
Uniprot	P34130
GeneID	4909;
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

## Application Details

Immunohistochemistry: 1:20 - 1:200

## Images



Immunohistochemical analysis of paraffin-embedded human prostate using #42517 at dilution of 1:100.

## Background

Target-derived survival factor for peripheral sensory sympathetic neurons.

## References

---

[1] "Cloning and characterization of cDNAs for murine macrophage inflammatory protein 2 and its human homologues." Tekamp-Olson P., Gallegos C., Bauer D., McClain J., Sherry B., Fabre M., van Deventer S., Cerami A.J. *Exp. Med.* 172:911-919(1990)

---

Note: This product is for in vitro research use only