

## Napsin-A Polyclonal Antibody

Catalog No: #42541

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Napsin-A Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous level of total Napsin-A polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant human Napsin-A protein
Target Name	Napsin-A
Other Names	Aspartyl protease 4 Napsin-1 TA01/TA02 NAPSA NAP1, NAPA
Accession No.	Swiss-Prot#: O96009
Uniprot	O96009
GeneID	9476;
Calculated MW	45kd
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

## Application Details

Western blotting:  1:500 - 1:1000

## Images



All lanes: Napsin-A antibody at 2ug/ml+mouse kidney tissue  
 Secondary  
 Goat polyclonal to Ribbit IgG at 1/10000 dilution  
 Predicted band size:45kDa  
 Observed band size:45kDa

## Background

May be involved in processing of pneumocyte surfactant precursors.

## References

---

[1]Napsins: new human aspartic proteinases. Distinction between two closely related genes.

---

Note: This product is for in vitro research use only