

Antigen KI-67 Polyclonal Antibody

Catalog No: #42544

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Antigen KI-67 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total Antigen KI-67 polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant human Antigen KI-67 protein
Target Name	Antigen KI-67
Other Names	MKI67
Accession No.	Swiss-Prot#: P46013
Uniprot	P46013
GeneID	4288;
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:20 - 1:200

Images



Immunohistochemical analysis of paraffin-embedded human endometrial cancer using #42544 at dilution of 1:100.



Immunohistochemical analysis of paraffin-embedded human bladder carcinoma using #42544 at dilution of 1:100.

Background

Thought to be required for maintaining cell proliferation.

References

[1] "The cell proliferation-associated antigen of antibody Ki-67: a very large, ubiquitous nuclear protein with numerous repeated elements, representing a new kind of cell cycle-maintaining proteins." Schlueter C., Duchrow M., Wohlenberg C., Be

Note: This product is for in vitro research use only