

Intestinal-type alkaline phosphatase Polyclonal Antibody

Catalog No: #42554

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Intestinal-type alkaline phosphatase Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	WB IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total Intestinal-type alkaline phosphatase polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant human Intestinal-type alkaline phosphatase protein
Target Name	Intestinal-type alkaline phosphatase
Other Names	IAP, Intestinal alkaline phosphatase, ALPI
Accession No.	Swiss-Prot#: P09923
Uniprot	P09923
GeneID	248;
Calculated MW	58kd
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

Application Details

Western blotting: □ 1:500 - 1:1000

Images



All lanes : Intestinal-type alkaline phosphatase antibody at 2ug/ml+293T whole cell lysate
 SecondaryGoat polyclonal to Rabbit IgG at 1/10000 dilution
 Predicted band size : 58 kDa
 Observed band size: 58 kDa

References

- [1] "Cloning and sequencing of human intestinal alkaline phosphatase cDNA."Berger J., Garattini E., Hua J.-C., Udenfriend S.Proc. Natl. Acad. Sci. U.S.A. 84:695-698(1987) [2] "Nucleotide and amino acid sequences of human intestinal alkaline phosphatase:

Note: This product is for in vitro research use only