

glsA1 Polyclonal Antibody

Catalog No: #42626

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	glsA1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	WB
Species Reactivity	E.coli
Specificity	The antibody detects endogenous level of total glsA1 polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant Escherichia coli Glutaminase 1 protein
Target Name	glsA1
Other Names	Corticotropin-releasing hormone receptor 1 CRHR1 CRFR, CRFR1, CRHR glsA1 ybaS b0485, JW0474
Accession No.	Swiss-Prot#: P77454
Uniprot	P77454
GeneID	946187;
Calculated MW	33kd
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

Application Details

Western blotting: □1:500 - 1:1000

Images



All lanes: Glutaminase 1 antibody at 2ug/ml+DH5α whole cell lysate
 secondary
 Goat polyclonal to rabbit at 1/10000 dilution
 predicted band size :33kDa
 observed band size :33kDa

References

[1]Sequence of minutes 4-25 of Escherichia coli.Chung E., Allen E., Araujo R., Aparicio A.M., Davis K., Duncan M., Federspiel N., Hyman R., Kalman S., Komp C., Kurdi O., Lew H., Lin D., Namath A., Oefner P., Roberts D., Schramm S., Davis R.W.Submitted (JA

Note: This product is for in vitro research use only