

CSRNP1 Polyclonal Antibody

Catalog No: #42663

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CSRNP1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen Affinity Purified
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total CSRNP1 polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant human Cysteine/serine-rich nuclear protein 1 protein(1-200aa)
Target Name	CSRNP1
Other Names	CSRNP-1, Axin-1 up-regulated gene 1 protein, Protein URAX1, TGF-beta-induced apoptosis protein 3, TAIP-3, CSRNP1, AXUD1, TAIP3
Accession No.	Swiss-Prot#: Q96S65
Uniprot	Q96S65
GeneID	64651;
Concentration	1.0mg/mL
Formulation	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:20 - 1:200

Images



Immunohistochemical analysis of paraffin-embedded human colon cancer using #42663 at dilution of 1:100.



Immunohistochemical analysis of paraffin-embedded human placenta using #42663 at dilution of 1:100.

Background

Binds to the consensus sequence 5'-AGAGTG-3' and has transcriptional activator activity (By similarity). May have a tumor-suppressor function. May play a role in apoptosis.

References

[1]The full-ORF clone resource of the German cDNA consortium."Bechtel S., Rosenfelder H., Duda A., Schmidt C.P., Ernst U., Wellenreuther R., Mehrle A., Schuster C., Bahr A., Bloecker H., Heubner D., Hoerlein A., Michel G., Wedler H., Koehrer K., Ottenwaelder B., Poustka A., Wiemann S., Schupp I.BMC Genomics 8:399-399(2007). [2]The status, quality, and expansion of the NIH full-length cDNA project: the Mammalian Gene Collection (MGC)."The MGC Project TeamGenome Res. 14:2121-2127(2004). [3]TGF-beta induced apoptosis protein 3 (TAIP-3)." Akiyama N., Kondoh S.Submitted (JUN-2001) .

Note: This product is for in vitro research use only