

Inactive phospholipase D5 Polyclonal Antibody

Catalog No: #42664

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Inactive phospholipase D5 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total Inactive phospholipase D5 polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant human Inactive phospholipase D5 protein
Target Name	Inactive phospholipase D5
Other Names	Inactive PLD 5, Inactive choline phosphatase 5, Inactive phosphatidylcholine-hydrolyzing phospholipase D5, PLDc, PLD5
Accession No.	Swiss-Prot#: Q8N7P1
Uniprot	Q8N7P1
GeneID	200150;
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:20 - 1:200

Images



Immunohistochemical analysis of paraffin-embedded human breast cancer using #42664 at dilution of 1:100.

References

[1]"Complete sequencing and characterization of 21,243 full-length human cDNAs."Ota T., Suzuki Y., Nishikawa T., Otsuki T., Sugiyama T., Irie R., Wakamatsu A., Hayashi K., Sato H., Nagai K., Kimura K., Makita H., Sekine M., Obayashi M., Nishi T., Shibahar

Note: This product is for in vitro research use only