

## Heat shock 70 kDa protein 6 Polyclonal Antibody

Catalog No: #42675

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Heat shock 70 kDa protein 6 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total Heat shock 70 kDa protein 6 polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant human Heat shock 70 kDa protein 6 protein
Target Name	Heat shock 70 kDa protein 6
Other Names	Heat shock 70 kDa protein B' HSPA6 HSP70B'
Accession No.	Swiss-Prot#: P17066
Uniprot	P17066
GeneID	3310;
Calculated MW	71kd
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

## Application Details

Western blotting: □1:500 - 1:1000

## Images



## Western blot

All lanes: Heat shock 70 kDa protein 6 antibody at 2ug/ml

Lane 1: MCF-7 whole cell lysate

Lane 2: A431 whole cell lysate

Lane 3: Jurkat whole cell lysate

secondary

Goat polyclonal to rabbit at 1/10000 dilution

predicted band size :71kDa

observed band size :71kDa

## References

[1]The human heat-shock protein family. Expression of a novel heat-inducible HSP70 (HSP70B') and isolation of its cDNA and genomic DNA. Leung T.K.C., Rajendran M.Y., Monfries C., Hall C., Lim L. Biochem. J. 267:125-132(1990) [2]The status, quality, and e

---

Note: This product is for in vitro research use only