

LMO4 Antibody

Catalog No: #42719

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	LMO4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total LMO4 protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human LMO4
Target Name	LMO4
Other Names	LIM domain only 4; LIM domain only protein 4; LIM domain transcription factor LMO4; LMO-4
Accession No.	Swiss-Prot#: P61968Gene ID: 8543
Uniprot	P61968
GeneID	8543;
Concentration	1.2mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human esophagus cancer tissue using #42719 at dilution 1/20,

Background

This gene encodes a cysteine-rich protein that contains two LIM domains but lacks a DNA-binding homeodomain. The encoded protein may play a role as a transcriptional regulator or as an oncogene.

Note: This product is for in vitro research use only