

TRIM69 Antibody

Catalog No: #42739

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TRIM69 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total TRIM69 protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human TRIM69
Target Name	TRIM69
Other Names	Trif; HSD34; RNF36; HSD-34
Accession No.	Swiss-Prot#: Q86WT6 Gene ID: 140691
Uniprot	Q86WT6
GeneID	140691;
Concentration	1.6mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human brain tissue using #42739 at dilution 1/30,

Background

This gene encodes a member of the RING-B-box-coiled-coil (RBCC) family and encodes a protein with an N-terminal RING finger motif, a PRY domain and a C-terminal SPRY domain. The mouse ortholog of this gene is specifically expressed in germ cells at the round spermatid stages during spermatogenesis and, when overexpressed, induces apoptosis. Alternatively spliced transcript variants encoding distinct isoforms have been described.?

Note: This product is for in vitro research use only