

SNX4 Antibody

Catalog No: #42765

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SNX4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total SNX4 protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human SNX4
Target Name	SNX4
Other Names	ATG24B
Accession No.	Swiss-Prot#: O95219 Gene ID: 8723
Uniprot	O95219
GeneID	8723;
Concentration	1.3mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human thyroid cancer tissue using #42765 at dilution 1/25,

Background

This gene encodes a member of the sorting nexin family. Members of this family contain a phox (PX) domain, which is a phosphoinositide binding domain, and are involved in intracellular trafficking. This protein associated with the long isoform of the leptin receptor and with receptor tyrosine kinases for platelet-derived growth factor, insulin, and epidermal growth factor in cell cultures, but its function is unknown. This protein may form oligomeric complexes with family members. Two transcript variants, one protein coding and the other non-protein coding, have been found for this

gene.

Note: This product is for in vitro research use only