

TFDP2 Antibody

Catalog No: #42777

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TFDP2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total TFDP2 protein.
Immunogen Type	protein
Immunogen Description	Full length fusion protein of human TFDP2
Target Name	TFDP2
Other Names	DP2
Accession No.	Swiss-Prot#: Q14188Gene ID: 7029
Uniprot	Q14188
GeneID	7029;
Concentration	1.8mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:30-1:150

Images



Immunohistochemical analysis of paraffin-embedded Human lung cancer tissue using #42777 at dilution 1/25,

Background

The gene is a member of the transcription factor DP family. The encoded protein forms heterodimers with the E2F transcription factors resulting in transcriptional activation of cell cycle regulated genes. Alternative splicing results in multiple transcript variants.

Note: This product is for in vitro research use only