

## TACSTD2 Antibody

Catalog No: #42791

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	TACSTD2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total TACSTD2 protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human TACSTD2
Target Name	TACSTD2
Other Names	EGP1; GP50; M1S1; EGP-1; TROP2; GA7331; GA733-1
Accession No.	Swiss-Prot#: P09758 Gene ID: 4070
Uniprot	P09758
GeneID	4070;
Concentration	2.1mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20°C

## Application Details

Immunohistochemistry: 1:25-1:100

## Images



Immunohistochemical analysis of paraffin-embedded Human breast cancer tissue using #42791 at dilution 1/35,

## Background

This intronless gene encodes a carcinoma-associated antigen. This antigen is a cell surface receptor that transduces calcium signals. Mutations of this gene have been associated with gelatinous drop-like corneal dystrophy.

---

Note: This product is for in vitro research use only