

TRPM2 Antibody

Catalog No: #42793

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TRPM2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total TRPM2 protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human TRPM2
Target Name	TRPM2
Other Names	KNP3; EREG1; TRPC7; LTRPC2; NUDT9H; NUDT9L1
Accession No.	Swiss-Prot#: O94759 Gene ID: 7226
Uniprot	O94759
GeneID	7226;
Concentration	1.4mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human liver cancer tissue using #42793 at dilution 1/25.



Immunohistochemical analysis of paraffin-embedded Human breast cancer tissue using #42793 at dilution 1/25.

Background

The protein encoded by this gene is a calcium-permeable cation channel that is regulated by free intracellular ADP-ribose. The encoded protein is activated by oxidative stress and confers susceptibility to cell death. Several alternatively spliced transcript variants of this gene have been described, but their full-length nature is not known.

Note: This product is for in vitro research use only