

TULP2 Antibody

Catalog No: #42797

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TULP2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total TULP2 protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human TULP2
Target Name	TULP2
Other Names	CT65; TUBL2
Accession No.	Swiss-Prot#: O00295Gene ID: 7288
Uniprot	O00295
GeneID	7288;
Concentration	1.2mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human liver cancer tissue using #42797 at dilution 1/20.



Immunohistochemical analysis of paraffin-embedded Human esophagus cancer tissue using #42797 at dilution 1/20.

Background

TULP2 is a member of a family of tubby-like genes (TULPs) that encode proteins of unknown function. Members of this family have been identified in plants, vertebrates, and invertebrates. The TULP proteins share a conserved C-terminal region of approximately 200 amino acid residues.?

Note: This product is for in vitro research use only