

TXNDC11 Antibody

Catalog No: #42799

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TXNDC11 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total TXNDC11 protein.
Immunogen Description	Fusion protein of human TXNDC11
Target Name	TXNDC11
Other Names	EFP1
Accession No.	Swiss-Prot#: Q6PKC3 Gene ID: 51061
Uniprot	Q6PKC3
GeneID	51061;
Concentration	2.7mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Liver cancer tissue using #42799 at dilution 1/100.

Background

May act as a redox regulator involved in DUOX proteins folding. The interaction with DUOX1 and DUOX2 suggest that it belongs to a multiprotein complex constituting the thyroid H₂O₂generating system. It is however not sufficient to assist DUOX1 and DUOX2 in H₂O₂?generation.

Note: This product is for in vitro research use only