

## UGGT2 Antibody

Catalog No: #42803

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	UGGT2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total UGGT2 protein.
Immunogen Type	protein
Immunogen Description	Full length fusion protein of human UGGT2
Target Name	UGGT2
Other Names	UGT2; HUGT2; UGCGL2
Accession No.	Swiss-Prot#: Q9NYU1 Gene ID: 55757
Uniprot	Q9NYU1
GeneID	55757;
Concentration	1.2mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20°C


## Application Details

Immunohistochemistry: 1:25-1:100

## Images



Immunohistochemical analysis of paraffin-embedded Human prostate cancer tissue using #42803 at dilution 1/20.



Immunohistochemical analysis of paraffin-embedded Human tonsil tissue using #42803 at dilution 1/20.

## Background

UDP-glucose:glycoprotein glucosyltransferase (UGT) is a soluble protein of the endoplasmic reticulum (ER) that selectively reglucosylates unfolded glycoproteins, thus providing quality control for protein transport out of the ER.

Note: This product is for in vitro research use only