

# Mouse IgG

Catalog No: #IG102



Package Size: #IG102-1 100mg #IG102-2 200mg #IG102-3 1g

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse IgG
Purification	Affinity chromatography or caprylic acid-ammonium sulfate,85%
Conjugates	None
Concentration	20 mg/ml
Storage	-20°C

## Product Description

Source: Normal, healthy, unimmunized mouse serum  
Appearance: Liquid  
Filtration: 0.2um  
Buffer :10mM Phosphate buffered saline  
pH7.40.  
Preservative: None  
Stability: Stable for 24 months  
Precautions: This product can be kept at 2-8°C for maximum 30 days. Avoid repeated freezing and thawing. If multiple usages are needed, it is recommended to be stored separately in smaller volumes.

Note: This product is for in vitro research use only