

Rabbit IgG

Catalog No: #IG401



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rabbit IgG
Purification	Affinity chromatography or caprylic acid-ammonium sulfate,95%
Conjugates	None
Concentration	mg/ml
Storage	-20°C

Product Description

Source: Normal, healthy, unimmunized rabbit serum

Appearance: Liquid

Filtration: 0.2um

Buffer :10mM Phosphate buffered salinePBSpH7.4

Preservative: None

Stability: Stable for 24 months

Precautions: This product can be kept at 2-8 Λ C for maximum 30 days. Avoid repeated freezing and thawing. If multiple usages are needed, it is recommended to be stored separately in smaller volumes.

Note: This product is for in vitro research use only