

Fc(human)IgG1Control(rec.)(non-lytic)

Catalog No: #AP79609



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Fc(human)IgG1Control(rec.)(non-lytic)
Brief Description	1545096460
Host Species	Human cell
Purification	>98%
Species Reactivity	Hu
Other Names	Fc region of human immunoglobulin G1
Calculated MW	34 KDa (SDS-PAGE)
Target Sequence	Mutant human IgG1 Fc region ,Pro100-Lys330 (Leu118Glu, Glu201Ala, Lys203Ala, Lys205Ala), containing the hinge sequence

Product Description

Endotoxin Content:

Note: This product is for in vitro research use only