

IL-12(human)(rec.)(His)

Catalog No: #AP79625



Package Size: #AP79625-1 50 ug #AP79625-2 10 ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IL-12(human)(rec.)(His)
Brief Description	1545096460
Host Species	Human cell
Purification	>95%
Species Reactivity	Hu
Other Names	CLMF; Cytotoxic lymphocyte maturation factor; NKSF; natural killer cell stimulatory factor; TSF; T-cell stimulating factor
Calculated MW	85kDa (SDS-PAGE)
Target Sequence	The recombinant IL-12 is comprised of p40 subunit (NP_002178.2) (ILe23-Ser328) and the p35 subunit (NP_000873.2) (Arg57-Ser253), these two subunit were linked by a (G3S)4 linker and with a C terminus 6xHis tag.

Product Description

Endotoxin Content:

Note: This product is for in vitro research use only