

# IL-3(human)(rec.)(His)

Catalog No: #AP79659



Package Size: #AP79659-1 50 ug #AP79659-2 10 ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | IL-3(human)(rec.)(His)  |
| Brief Description  | 1545096460  |
| Host Species       | Human cell  |
| Purification       | >95%  |
| Species Reactivity | Hu  |
| Other Names        | IL-3; Interleukin-3   |
| Calculated MW      | 20kDa (SDS-PAGE)  |
| Target Sequence    | Human IL-3 (AAH69472.1) (Ala20-Phe152), with a C-terminal 6-His tag |

## Product Description

Endotoxin Content:

Note: This product is for in vitro research use only