

VPS25 Antibody

Catalog No: #42837

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	VPS25 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total VPS25 protein.
Immunogen Type	protein
Immunogen Description	Full length fusion protein of human VPS25
Target Name	VPS25
Other Names	DERP9; EAP20; FAP20
Accession No.	Swiss-Prot#: Q9BRG1 Gene ID: 84313
Uniprot	Q9BRG1
GeneID	84313;
Calculated MW	21kd
Concentration	1.3mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:200-1:1000

Immunohistochemistry: 1:25-1:100

Images



Gel: 12%SDS-PAGE
Lysate: 40 µg
Lane: Human fetal liver tissue
Primary antibody: 1/200 dilution
Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution
Exposure time: 30 seconds

Background

This gene encodes a protein that is a subunit of the endosomal sorting complex required for transport II (ESCRT-II). This protein complex functions in

sorting of ubiquitinated membrane proteins during endocytosis. A pseudogene of this gene is present on chromosome 1.?

Note: This product is for in vitro research use only