

IL-4(mouse)(rec.)(His)

Catalog No: #80324



Package Size: #80324-1 10 ug #80324-2 50 ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	IL-4(mouse)(rec.)(His)
Brief Description	1545096460
Host Species	E. coli
Purification	>98%
Species Reactivity	Ms
Other Names	Interleukin-4, B Cell Stimulatory Factor 1,BSF-1, BCGF-1, BCGF1
Calculated MW	15KDa
Target Sequence	Mouse IL-4(NP_067258.1)(His21-Ser140) is fused at the C-terminus to a His-tag.

Product Description

Endotoxin Content:

Note: This product is for in vitro research use only