

TIM3 Antibody

Catalog No: #P1003

Package Size: #P1003-1 0.1ml #P1003-2 1ml

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	TIM3 Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	IHC003
Applications	IHC
Species Reactivity	Human
Formulation	Tris Buffer, pH 7.3 - 7.7, with 1% BSA and
Storage	Store at 2-8C. Do not freeze.

Application Details

Recommended

working dilution range:1:100 - 1:200

Images



TIM3 Antibody (P1003) on Kidney - 40X

Product Description

T-cell Immunoglobulin and Mucin-Domain-Containing Molecule-3 (TIM3) is present on T-helper type 1 lymphocytes and other immune cells including dendritic cells and natural killer cells. TIM3 is overexpressed in CD4+ tumor-infiltrating lymphocytes, including those with non-small cell lung cancer associated with poor prognoses. TIM3 has recently emerged as a potential target for cancer immunotherapy.

References

- Monney L, et al. Nature. 2002; 415:536-41.
- Snchez-Fueyo A, et al. Nat Immunol. 2003; 4:1093-101.
- Anderson AC. Cancer Immun Res. 2014; 2:393C8.
- Gao X, et al. PLoS One. 2012; 7:e30676.

Note: This product is for in vitro research use only