

## TAG-72 Antibody

Catalog No: #P1023

Package Size: #P1023-1 0.1ml #P1023-2 1ml

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	TAG-72 Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	IHC072
Applications	IHC
Species Reactivity	Human
Formulation	Tris Buffer, pH 7.3 - 7.7, with 1% BSA and
Storage	Store at 2-8C. Do not freeze.

## Application Details

Recommended

working dilution range:1:100 - 1:200

## Images



TAG-72 Antibody (P1023) on Colon

## Product Description

Tumor-Associated Glycoprotein 72 (TAG-72) is a glycoprotein found on the surface of many cancer pathologies. Anti-TAG-72 can be useful for detecting some adenocarcinomas and non-neoplastic tissues. This diagnostic grade TAG-72 IVD antibody is useful for identifying adenocarcinomas in positive staining, but in mesotheliomas no staining is observed.

## References

- Thor A, et al. *Cancer Res.* 1986; 46:3118-24.  
Osteen KG, et al. In *J Gynecol Pathol.* 1992; 11:216-20.  
Johnston WW, et al. *Hum Pathol.* 1986; 17:501-13.  
Lundy J, et al. *Ann Surg.* 1986; 203:399-402.  
Kline TS, et al. *Cancer.* 1989; 63:2253-6.

Chieng DC, et al. Hum Pathol. 2003; 34:1016-21.  
Ord?ez NG. Am J Surg Pathol. 1998; 22:1203-14.  
Ord?ez NG. Am J Surg Pathol. 2003; 27:1031-51.

---

Note: This product is for in vitro research use only