

## Vimetin Antibody

Catalog No: #P1055

Package Size: #P1055-1 0.1ml #P1055-2 1ml

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	Vimetin Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	IHC684
Applications	IHC
Species Reactivity	Human
Formulation	Tris Buffer, pH 7.3 - 7.7, with 1% BSA and
Storage	Store at 2-8C. Do not freeze.

## Application Details

Recommended

working dilution range:1:100 - 1:200

## Images



Vimentin [P1055] on Colon

## Product Description

Vimentin is a component of intermediate filament in mesenchymal cells, such as endothelial cells, fibroblasts, lymphocytes, and melanocytes. Anti-Vimentin is useful for assessing whether tissue samples have been processed and preserved properly. A panel of Anti-Vimentin and Anti-Keratin is useful for differentiating melanomas from large cell lymphomas and undifferentiated carcinomas. This diagnostic grade Vimentin IVD antibody stains melanomas and schwannomas, as well as Endometrial endometrioid adenocarcinomas.

## References

- Dabbs DJ, et al. Hum Pathol. 1996; 27:172-7.  
 Dabbs DJ, et al. Am J Surg Pathol. 1986; 10:568-76.  
 Yaziji H, et al. Int J Gynecol Pathol. 2001; 20:64-78.

---

Note: This product is for in vitro research use only