

# Recombinant Human Myeloblastin(PRTN3),Escherichia Coli

Catalog No: #AP72188

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP72188-1 20ug #AP72188-2 100ug #AP72188-3 1mg

## Description

Product Name	Recombinant Human Myeloblastin(PRTN3),Escherichia Coli
Host Species	E.coli
Purification	>85% (SDS-PAGE)
Other Names	MBN
Accession No.	Uniprot ID:P24158 Gene Names:PRTN3
Uniprot	P24158
GeneID	5657;
Tag Info	N-termina His-tagged
Target Sequence	YMASLQMRGN PGSHFCGGTL IHPSFVLTAH HCLRDIPQRL VNVVLGAHNV RTQEPTQQHF SVAQVFLNNY DAENKLNVDL LIQLSSPANL SASVATVQLP QQDQPVPHGT QCLAMGWGRV GAHDPPAQVL QELNVTVVTF FCRPHNICTF VPRRKAGICF GDSGGP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

## Product Description

Note: Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Storage Buffer: Tris-based buffer, 50% glycerol

Organism: Homo sapiens (Human)

Expression Region: 28-248aa

Sequence Info: Full Length

Note: This product is for in vitro research use only