

Recombinant Human HSP90AA1

Catalog No: #GP10006

Package Size: #GP10006-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human HSP90AA1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250 amino acids of human heat shock protein 90kDa alpha (cytosolic), class A member 1
Target Name	heat shock protein 90kDa alpha (cytosolic), class A member 1
Other Names	EL52; HSPN; LAP2; HSP86; HSPC1; HSPCA; Hsp89; Hsp90; LAP-2; HSP89A; HSP90A; HSP90N; HSPCAL1; HSPCAL4
Accession No.	Swissprot:P07900Gene Accession:BC121062
Uniprot	P07900
GeneID	3320;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is an inducible molecular chaperone that functions as a homodimer. The encoded protein aids in the proper folding of specific target proteins by use of an ATPase activity that is modulated by co-chaperones. Two transcript variants encoding different isoforms have been found for this gene

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only