

Recombinant Human TFF3

Catalog No: #GP10009



Package Size: #GP10009-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human TFF3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 22-80 amino acids of human TFF3
Target Name	trefoil factor 3 (intestinal)
Other Names	ITF; P1B; TFI
Accession No.	Swissprot:Q07654Gene Accession:BC017859
Uniprot	Q07654
GeneID	7033;
Storage	-20~-80°C, pH 7.6 PBS

Background

Members of the trefoil family are characterized by having at least one copy of the trefoil motif, a 40-amino acid domain that contains three conserved disulfides. They are stable secretory proteins expressed in gastrointestinal mucosa. Their functions are not defined, but they may protect the mucosa from insults, stabilize the mucus layer and affect healing of the epithelium. This gene is expressed in goblet cells of the intestines and colon.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only