

Recombinant Human UPK3B

Catalog No: #GP10016



Package Size: #GP10016-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human UPK3B
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	Uroplakin-3b
Other Names	FLJ32198, MGC10902, p35, UPIIIb
Accession No.	Swissprot:Q9BT76Gene Accession:BC004304
Uniprot	Q9BT76
GeneID	105375355;
Storage	-20~-80°C, pH 7.6 PBS

Background

UPK3B is a minor component of the apical plaques of mammalian urothelium that binds and dimerizes with uroplakin-1b (UPK1B, MIM 602380), one of the major conserved urothelium membrane proteins. The other major conserved integral membrane proteins of urothelial plaques are UPK1A (MIM 611557), UPK2 (MIM 611558), and UPK3A (MIM 611559)

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only