

# Recombinant Human ADAP1

Catalog No: #GP10039



Package Size: #GP10039-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ADAP1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 3-376 amino acids of human Arf-GAP with dual PH domain-containing protein 1
Target Name	Arf-GAP with dual PH domain-containing protein 1
Other Names	GCS1L, CENTA1, p42IP4
Accession No.	Swissprot:O75689Gene Accession:BC033747
Uniprot	O75689
GeneID	11033;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Arf-GAP with dual PH domain-containing protein 1 is a protein that in humans is encoded by the ADAP1 gene. GTPase-activating protein for the ADP ribosylation factor family. Binds phosphatidylinositol 3,4,5-trisphosphate (PtdInsP3) and inositol 1,3,4,5-tetrakisphosphate (InsP4).

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only