

# Recombinant Human JUN

Catalog No: #GP10054



Package Size: #GP10054-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human JUN
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250 amino acids of human jun proto-oncogene
Target Name	jun proto-oncogene
Other Names	AP1; AP-1; c-Jun
Accession No.	Swissprot:P05412Gene Accession:BC006175
Uniprot	P05412
GeneID	3725;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene is the putative transforming gene of avian sarcoma virus 17. It encodes a protein which is highly similar to the viral protein, and which interacts directly with specific target DNA sequences to regulate gene expression. This gene is intronless and is mapped to 1p32-p31, a chromosomal region involved in both translocations and deletions in human malignancies.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only