

DPH7 Antibody

Catalog No: #42855

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DPH7 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total DPH7 protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human DPH7
Target Name	DPH7
Other Names	RRT2; WDR85; C9orf112
Accession No.	Swiss-Prot#: Q9BTV6 Gene ID: 92715
Uniprot	Q9BTV6
GeneID	92715;
Concentration	2.2mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human lung cancer tissue using #42855 at dilution 1/40.



Immunohistochemical analysis of paraffin-embedded Human thyroid cancer tissue using #42855 at dilution 1/40.

Background

WDR85 is a WD repeat-containing protein that plays a role in the first step of diphthamide biosynthesis (Carette et al., 2009 [PubMed 19965467]).

Note: This product is for in vitro research use only