

Recombinant Human COX16

Catalog No: #GP10090



Package Size: #GP10090-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human COX16
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	COX16 cytochrome c oxidase assembly homolog (S. cerevisiae)
Other Names	HSPC203, C14orf112
Accession No.	Swissprot:Q9P0S2Gene Accession:BC001702
Uniprot	Q9P0S2
GeneID	51241;
Storage	-20~-80°C, pH 7.6 PBS

Background

This protein is encoded by COX16, corresponding to the previously uncharacterized open reading frame YJL003w of the yeast genome. COX16 was identified in studies of COX-deficient mutants previously assigned to complementation group G22 of a collection of yeast pet mutants. To determine its location, Cox16p was tagged with a Myc epitope at the C terminus.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only