

Recombinant Human COX19

Catalog No: #GP10096

Package Size: #GP10096-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human COX19
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	COX19 cytochrome c oxidase assembly factor
Other Names	MGC104475
Accession No.	Swissprot:Q49B96Gene Accession:BC110420
Uniprot	Q49B96
GeneID	90639;
Storage	-20~-80°C, pH 7.6 PBS

Background

COX19 encodes a cytochrome c oxidase (COX)-assembly protein. The *S. cerevisiae* Cox19 protein may play a role in metal transport to the mitochondrial intermembrane space and assembly of complex IV of the mitochondrial respiratory chain.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only