

Recombinant Human IL20RA

Catalog No: #GP10133



Package Size: #GP10133-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human IL20RA
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 30-250 amino acids of human interleukin 20 receptor, alpha
Target Name	interleukin 20 receptor, alpha
Other Names	IL-20R1; ZCYTOR7
Accession No.	Swissprot:Q9UHF4Gene Accession:BC113574
Uniprot	Q9UHF4
GenID	53832;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is a receptor for interleukin 20 (IL20), a cytokine that may be involved in epidermal function. The receptor of IL20 is a heterodimeric receptor complex consisting of this protein and interleukin 20 receptor beta (IL20B). This gene and IL20B are highly expressed in skin. The expression of both genes is found to be upregulated in Psoriasis.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only