

Recombinant Human ABCB5

Catalog No: #GP10143



Package Size: #GP10143-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ABCB5
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1007-1257 amino acids of human ABCB5
Target Name	ATP binding cassette subfamily B member 5
Other Names	ABCB5beta; EST422562; ABCB5alpha
Accession No.	Swissprot:Q2M3G0Gene Accession:BC104894
Uniprot	Q2M3G0
GeneID	340273;
Storage	-20~-80°C, pH 7.6 PBS

Background

ABCB5 belongs to the ATP-binding cassette (ABC) transporter superfamily of integral membrane proteins. These proteins participate in ATP-dependent transmembrane transport of structurally diverse molecules ranging from small ions, sugars, and peptides to more complex organic molecules (Chen et al., 2005 [PubMed 15760339]).

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only