

Recombinant Human AMPH

Catalog No: #GP10155



Package Size: #GP10155-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human AMPH
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250 amino acids of human amphiphysin
Target Name	amphiphysin
Other Names	AMPH1
Accession No.	Swissprot:P49418Gene Accession:BC034376
Uniprot	P49418
GeneID	273;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a protein associated with the cytoplasmic surface of synaptic vesicles. A subset of patients with stiff-man syndrome who were also affected by breast cancer are positive for autoantibodies against this protein. Alternate splicing of this gene results in two transcript variants encoding different isoforms. Additional splice variants have been described, but their full length sequences have not been determined. A pseudogene of this gene is found on chromosome 11.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only