

Recombinant Human ANGPTL1

Catalog No: #GP10156



Package Size: #GP10156-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ANGPTL1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 143-491 amino acids of human angiopoietin-like 1
Target Name	angiopoietin-like 1
Other Names	ANG3; ARP1; AngY; ANGPT3; UNQ162; dJ595C2.2
Accession No.	Swissprot:O95841Gene Accession:BC050640
Uniprot	O95841
GeneID	9068;
Storage	-20~-80°C, pH 7.6 PBS

Background

Angiopoietins are members of the vascular endothelial growth factor family and the only known growth factors largely specific for vascular endothelium. Angiopoietin-1, angiopoietin-2, and angiopoietin-4 participate in the formation of blood vessels. The protein encoded by this gene is another member of the angiopoietin family that is widely expressed in adult tissues with mRNA levels highest in highly vascularized tissues. This protein was found to be a secretory protein that does not act as an endothelial cell mitogen in vitro.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only