

XAB2 Antibody

Catalog No: #42863

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	XAB2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total XAB2 protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human XAB2
Target Name	XAB2
Other Names	HCNP; HCRN; SYF1; NTC90
Accession No.	Swiss-Prot#: Q9HCS7 Gene ID: 56949
Uniprot	Q9HCS7
GeneID	56949;
Concentration	0.5mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human lung cancer tissue using #42863 at dilution 1/20.



Immunohistochemical analysis of paraffin-embedded Human thyroid cancer tissue using #42863 at dilution 1/20.

Background

Involved in transcription-coupled repair (TCR), transcription and pre-mRNA splicing.

Note: This product is for in vitro research use only