

Recombinant Human ANGPTL3

Catalog No: #GP10159



Package Size: #GP10159-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ANGPTL3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 161-460 amino acids of human ANGPTL3
Target Name	angiopoietin-like 3
Other Names	ANL3; ANG-5; FHBL2; ANGPT5
Accession No.	Swissprot:Q9Y5C1 Gene Accession:BC058287
Uniprot	Q9Y5C1
GeneID	27329;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of a family of secreted proteins that function in angiogenesis. The encoded protein, which is expressed predominantly in the liver, is further processed into an N-terminal coiled-coil domain-containing chain and a C-terminal fibrinogen chain. The N-terminal chain is important for lipid metabolism, while the C-terminal chain may be involved in angiogenesis. Mutations in this gene cause familial hypobetalipoproteinemia type 2.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only