

Recombinant Human ANGPT2

Catalog No: #GP10162



Package Size: #GP10162-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ANGPT2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 300 amino acids of human angiotensin 2
Target Name	angiotensin 2
Other Names	ANG2; AGPT2
Accession No.	Swissprot:O15123Gene Accession:BC126200
Uniprot	O15123
GeneID	285;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is an antagonist of angiotensin 1 (ANGPT1) and endothelial TEK tyrosine kinase (TIE-2, TEK). The encoded protein disrupts the vascular remodeling ability of ANGPT1 and may induce endothelial cell apoptosis. Three transcript variants encoding three different isoforms have been found for this gene.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only