

# Recombinant Human ANXA1

Catalog No: #GP10163



Package Size: #GP10163-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ANXA1
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	Annexin A1
Other Names	ANX1; LPC1
Accession No.	Swissprot:P04083Gene Accession:BC001275
Uniprot	P04083
GeneID	301;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Annexin I belongs to a family of Ca(2+)-dependent phospholipid binding proteins which have a molecular weight of approximately 35,000 to 40,000 and are preferentially located on the cytosolic face of the plasma membrane. Annexin I protein has an apparent relative molecular mass of 40 kDa, with phospholipase A2 inhibitory activity. Since phospholipase A2 is required for the biosynthesis of the potent mediators of inflammation, prostaglandins and leukotrienes, annexin I may have potential anti-inflammatory activity.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only