

# Recombinant Human ANXA3

Catalog No: #GP10164



Package Size: #GP10164-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ANXA3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-322 amino acids of human annexin A3
Target Name	annexin A3
Other Names	ANX3
Accession No.	Swissprot:P12429Gene Accession:BC000871
Uniprot	P12429
GeneID	306;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a member of the annexin family. Members of this calcium-dependent phospholipid-binding protein family play a role in the regulation of cellular growth and in signal transduction pathways. This protein functions in the inhibition of phospholipase A2 and cleavage of inositol 1,2-cyclic phosphate to form inositol 1-phosphate. This protein may also play a role in anti-coagulation.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only