

# Recombinant Human ANXA8

Catalog No: #GP10167



Package Size: #GP10167-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ANXA8
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 130-200 amino acids of human annexin A8
Target Name	annexin A8
Other Names	ANX8
Accession No.	Swissprot:P13928Gene Accession:BC073755
Uniprot	P13928
GeneID	653145;728113;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a member of the annexin family of evolutionarily conserved Ca<sup>2+</sup> and phospholipid binding proteins. The encoded protein may function as an anticoagulant that indirectly inhibits the thromboplastin-specific complex. Overexpression of this gene has been associated with acute myelocytic leukemia. A highly similar duplicated copy of this gene is found in close proximity on the long arm of chromosome 10.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only